

Optimus IMF Screw



OSTEONIC
The Medical Explorer



Introduction

Optimus IMF Screw system is the most optimized for clinicians

Intermaxillary fixation (IMF) Screw gives surgeons a fast and effective alternative to arch bars.

The hole just under head of screw is designed for convenient use of standard ligating wire.

Securing easy and fast reach of IMF screw, it is installed at both Maxillofacial and Mandible sets.

D-Head(Without hole)



No.	Item	D(mm)	L(mm)
01	D21606	1.6	6
02	D21608		8
03	D21610		10
04	D21612		12
05	D21806	1.8	6
06	D21808		8
07	D21810		10
08	D21812		12
09	D22006	2.0	6
10	D22008		8
11	D22010		10
12	D22012		12

D-Head(With hole)



No.	Item	D(mm)	L(mm)
13	D31606	1.6	6
14	D31608		8
15	D31610		10
16	D31612		12
17	D31806	1.8	6
18	D31808		8
19	D31810		10
20	D31812		12
21	D32006	2.0	6
22	D32008		8
23	D32010		10
24	D32012		12

Instrument

Item	Item Code	Description
	IDS-01	Driver Shaft 70mm for D-head Type
	IDS-02	Driver Shaft 40mm for D-head Type
	IDS-03	Driver Shaft 23mm for D-head Type

- Simplified intermaxillary fixation eliminates arch bars
- Reduces treatment time by drill free screw design
- Innovative screw design allows for low insertion torque.

H-Head



No.	Item	D(mm)	L(mm)
32	D42006B	2.0	6
33	D42007B		7
34	D42008B		8
35	D42009B		9
36	D42010B		10
37	D42011B		11
38	D42012B		12

Instrument

Item	Item Code	Description
	ODS-01	Cruciform Driver Shaft 90 H-head Type
	IMF-L04	Driver Shaft 70mm for H-head Type

CE
1639



OSTEONIC Co., Ltd.

1206Ho, 38, Digital-ro 29-gil, Guro-gu, Seoul, Korea P.O. #08381
TEL: +822 6902 8400 FAX: +822 6902 8401 www.osteonic.com

IFS - E - 2005